Veterinary Endoscope for Insemination & Hysteroscopy





A Revolution in Artificial Insemination

Endoscopy is the best method for diagnosing and treating conditions of the genitourinary tract in bitches. The fine diameter and extended length of the KARL STORZ TCI endoscope allows easy access to the cervix and uterus in bitches of all sizes. An integrated working channel enables the insertion of instruments under direct visualization for procedures such as insemination, biopsy, aspiration, flushing and removal of foreign material.



Prof. Alain FONTBONNE, DVM, PhD, Dip.ECAR, Head of the Centre d'Etude en Reproduction des Carnivores (CERCA) at Alfort National Veterinary College in Paris, France, performs transcervical insemination in a bitch with assistance from Paulo Borges, DVM and Cindy Maenhoudt, DVM, DACT, DECAR. For a listing of available seminars, visit karlstorz.com or contact us at vetglobal@karlstorz.com

Compared to surgical insemination, endoscopic transcervical insemination (TCI) is less painful, less invasive and no post-operative care is needed. TCI is typically performed without sedation, making it easier than other methods.

The KARL STORZ TCI endoscope is also suitable for transcervical insemination in small ruminants and large non-domestic felids. For information about laparoscopic artificial insemination in small ruminants, see SELLSHEET VET 4: *Laparoscopic Artificial Insemination in Sheep & Goats.*

Applications:

- Visual evaluation of the vagina, cervix and uterus
- Transcervical Insemination (TCI) possibility of several inseminations during the fertile period
- Diagnosis and possible treatment of vaginal defects (septa, strictures, adhesions)
- Visualization of the vaginal wall and cervix during TCI, addressing possible problems
 (consistive) to improve TCI
- (cervicitis) to improve TCI outcomes
- Biopsy of masses in the vagina, cervix and uterus
- Flushing of the uterus to collect samples
- Diagnosis of a retained fetus or placenta (and possible removal)
- Assessment of obstructive dystocia




Insertion of the endoscope into the vagina. The use of a finger over the end of the endoscope during insertion may help to prevent inadvertently introducing the telescope into the urethral opening.









- Fig.1: Catheterization of the cervix for artificial insemination or diagnostic procedures, including uterine biopsies
- Fig.2: Hysteroscopy: presence of petechiae on the endometrium (possible sign of endometritis)
- Fig.3: Vaginal wall during oestrus, showing the crenulations
- Fig.4: Visualization of an amniotic sac in the cranial vagina of a bitch at stage 1 of parturition


Hysteroscope/TCI Endoscope

Special features:

- Atraumatic distal sheath with rounded tip
- Minimal shaft diameter, extended length facilitates cervical catheterization in small and large patients
- The large working channel allows the use of rigid instruments and catheters up to 5 Fr.
- 2 right angled lateral irrigation connections
- The rigid angled eyepiece makes working more comfortable and ensures easy introduction of catheters and instruments without bending or kinking

Specifications:

Distal tip:	8 Fr.
Sheath:	9.5 Fr., conical, 1 step 9.5-13.5 Fr.
Working channel:	6 Fr., for use with instruments up to 5 Fr.
Telescope:	KARL STORZ Fiberoptic system, direction of view 6
Length:	43 cm
Eyepiece:	angled, rigid





63002LA **Hys**

Hysteroscope/TCI Endoscope, length 43 cm, autoclavable



The following accessories are included:

27001G	

Instrument Port with Sealing System and Quick Release Lock
Seal, for Instrument Port 27001G, package of 10
Cleaning Adaptor, for port 27001G

Optional Instruments

	27425F	Forceps for grasping stone fragments, rigid, double action jaws, 5 Fr., length 60 cm
\rightarrow	27425P	Forceps for grasping larger stones and fragments, rigid, double action jaws, 5 Fr., length 60 cm
2	27425Z	Biopsy Forceps, rigid, double action jaws, 5 Fr., length 60 cm
	27651UB	Cleaning Brush, round, flexible, outer diameter 2.2 mm, for working channel diameter 0.8-2 mm, length 60 cm
	63002MK	TCI Catheter, 5 Fr., length 70 cm

Wire Tray for Cleaning, Sterilization and Storage




39501X Wire Tray for Cleaning, Sterilization and Storage, with integrated cleaning adaptor for washer-disinfector, with lid, spare parts basket 39501XS and silicone telescope holders, external dimensions (w x d x h): 644 x 150 x 80 mm



Check out our schedule of upcoming hands-on training courses at http://go.karlstorz.com/eventsVET

Notes



Notes

	_													



Shaping the Future of Endoscopy with you



KARL STORZ SE & Co. KG Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany Postbox 230, 78503 Tuttlingen/Germany Phone: +49 7461 708-0 Fax: +49 7461 708-105 E-Mail: info@karlstorz.com

www.karlstorz.com